

SYSTEM FLOW CHART

**MAQUET**  
GETINGE GROUP

ANESTHESIA  
PATIENT CONNECTION



# PATIENT CONNECTION ANESTHESIA

## STARTER KIT CONSUMABLES

One selected kit may be ordered with every FLOW-i

### 22 mm patient tubes

**Starter kit Anesthesia, adult, 22 mm**  
66 96 810


- Including:
  - 1 pc of Set of breathing circuit, 22 mm, disposable (20 pcs: 66 92 205).
  - 1 pc of Servo Duo Guard viral/bacterial filter (60 pcs: 66 71 775).
  - 1 pc of Anesthesia mask, disposable, Size 4, medium adult (20 pcs: 66 92 242).

### 10 mm patient tubes

**Starter kit anesthesia, infant, 10 mm mm**  
66 96 712

- Including:
  - 1 pc of Set of breathing circuit, 10 mm, disposable (20 pcs: 66 92 210).
  - 1 pc of Servo Duo Guard viral/bacterial filter (60 pcs: 66 71 775).
  - 1 pc of Anesthesia mask, reusable, Size 0, small infant (10 pcs: 66 92 225).

Note! Not all products listed here are available in all countries. For details contact your sales representative.



## HME Disposable

### HEAT AND MOISTURE EXCHANGE (HME)



Regular volume	Regular volume	Without HME	Regular volume	Regular volume	Small volume
<p><b>Servo Humidifier 162, 100 pcs/pkg</b> 64 19 357</p>	<p><b>Servo Humidifier 163, 100 pcs/pkg</b> With flex tube 64 19 365</p>	<p><b>Servo Filter Humidifier 172, 100 pcs/pkg</b> 64 19 373</p>	<p><b>Servo Filter Humidifier 173, 100 pcs/pkg</b> With flex tube 64 19 381</p>	<p><b>Servo Humidifier 161, 100 pcs/pkg</b> 64 19 340</p>	

## PATIENT CONNECTION


Reusable			Disposable		
22 mm patient tubes	15 mm patient tubes	10 mm patient tubes	22 mm patient tubes	15 mm patient tubes	10 mm patient tubes
<p><b>Breathing circuit, 22 mm, Hytrel (1.5 m)</b> 66 92 180</p> <ul style="list-style-type: none"> <li>• Including:                             <ul style="list-style-type: none"> <li>- 2 x Patient tube 150 cm (22 mm conn).</li> <li>- Y-piece 22 mm conn.</li> <li>- Angled adapter 22/15 mm conn. with luer-lock port.</li> </ul> </li> </ul>	<p><b>Breathing circuit, 15 mm, Hytrel (1.5 m)</b> 66 92 182</p> <ul style="list-style-type: none"> <li>• Including:                             <ul style="list-style-type: none"> <li>- 2 x Patient tube 150 cm (22/15 mm conn).</li> <li>- Y-piece 22/15 mm conn.</li> <li>- Angled adapter 22/15 mm conn. with luer-lock port.</li> </ul> </li> </ul>	<p><b>Breathing circuit, 10 mm, Hytrel (1.5 m)</b> 66 92 185</p> <ul style="list-style-type: none"> <li>• Including:                             <ul style="list-style-type: none"> <li>- 2 x Patient tube 150 cm (22/10 mm conn).</li> <li>- Angled Y-piece with luer-lock port, 15/10 mm conn.</li> </ul> </li> </ul>	<p><b>Set of breathing circuit, 22 mm, disposable (1.5 m), 20 pcs/pkg</b> 66 92 205</p> <ul style="list-style-type: none"> <li>• Including:                             <ul style="list-style-type: none"> <li>- Patient tubes (smooth-bore), 22 mm diameter, 150 cm.</li> <li>- Manual bag tube (smooth-bore), 22 mm, 100 cm.</li> <li>- Latex-free manual bag, 2 L.</li> <li>- Y-piece.</li> <li>- Gas sampling line, 2.0 m.</li> <li>- Elbow connector with gas sampling port.</li> </ul> </li> </ul>	<p><b>Set of breathing circuit, 15 mm, disposable (1.2 m), 20 pcs/pkg</b> 66 92 207</p> <ul style="list-style-type: none"> <li>• Including:                             <ul style="list-style-type: none"> <li>- Patient tubes (smooth-bore), 15 mm diameter with 22 mm connectors, 120 cm.</li> <li>- Manual bag tube (smooth-bore), 15 mm diameter with 22 mm connector, 120 cm.</li> <li>- Latex-free manual bag, 1 L.</li> <li>- Y-piece.</li> <li>- Gas sampling line, 2.0 m.</li> <li>- Elbow connector with gas sampling port.</li> </ul> </li> </ul>	<p><b>Set of breathing circuit, 10 mm, disposable (1.2 m), 20 pcs/pkg</b> 66 92 210</p> <ul style="list-style-type: none"> <li>• Including:                             <ul style="list-style-type: none"> <li>- Patient tubes (smooth-bore), 10 mm diameter with 22 mm connectors, 120 cm.</li> <li>- Manual bag tube (smooth-bore), 10 mm diameter with 22 mm connector, 120 cm.</li> <li>- Latex-free manual bag, 0.5 L.</li> <li>- Y-piece.</li> <li>- Gas sampling line, 2.0 m.</li> <li>- Elbow connector with gas sampling port.</li> </ul> </li> </ul>
<p><b>Reusable</b></p> <p><b>Tube connector 22-15 mm, angled, 5 pcs/pkg</b> 90 37 730</p>			<p><b>Reusable</b></p> <p><b>Tube joint/adapter 22 mm</b> 64 25 581</p> <ul style="list-style-type: none"> <li>• 22 mm out. diam., 15 mm int. diam.</li> </ul>		

**ANESTHESIA MASK**

**Reusable**


	<b>Anesthesia mask, reusable, 10 pcs/pkg</b>	
	Size 0, small infant	66 92 225
	Size 1, medium infant	66 92 227
	Size 2, large infant	66 92 230
	Size 3, small adult	66 92 232
	Size 4, medium adult	66 92 235
	Size 5, large adult	66 92 237

**Disposable**

	<b>Anesthesia mask, disposable, 20 pcs/pkg</b>	
	Size 3, small adult	66 92 240
	Size 4, medium adult	66 92 242
	Size 5, large adult	66 92 245
	Size 6, XL adult	66 92 247

**VIRAL/BACTERIAL FILTER**


**Disposable**

	<b>Servo DuoGuard viral/bacterial filter, 60 pcs/pkg</b>	
		66 71 775


**MANUAL VENTILATION**

**Reusable**


**22 mm patient tubes**

	<b>Breathing bag set 22 mm, Hytrel</b>	
	66 92 187	
	• Including:	
	- Patient tube 110 cm (22 mm conn).	
	- Breathing bag 2 L (reusable).	
	- Adapter 22 mm conn.	

**15 mm patient tubes**


	<b>Breathing bag set 15 mm, Hytrel</b>	
	66 92 192	
	• Including:	
	- Patient tube 110 cm (22 mm conn).	
	- Breathing bag 1 L (reusable).	
	- Adapter 22 mm conn.	

**10 mm patient tubes**


	<b>Breathing bag set 10 mm, Hytrel</b>	
	66 92 190	
	• Including:	
	- Patient tube 110 cm (22 mm conn).	
	- Breathing bag 0.5 L (reusable).	
	- Adapter 22 mm conn.	

**Alternative breathing bag**

**Reusable**


	<b>Breathing bag, reusable</b>	
	0.5 L	66 92 195
	1 L	66 92 197
	2 L	66 92 200
	3 L	66 92 202

**Disposable**

	<b>Breathing bag, disposable, 20 pcs/pkg</b>	
	0.5 L	66 92 215
	1 L	66 92 217
	2 L	66 92 220
	3 L	66 92 222


**O<sub>2</sub> FLOWMETER / SUCTION EQUIPMENT**

**Extra filter Disposable**

	<b>Hydrofobic bacterial/viral filter</b>	
	57 500 630 *	
	• For Auxiliary O <sub>2</sub> and suction module.	
	• Disposable.	
* Order from Maquet KG, Rastatt, Business Unit MEDAP, phone +49-7222-932-0.		



## GAS ANALYZE

### Disposable




## GAS EVACUATION

### Disposable




## CO<sub>2</sub> ABSORBER


### Disposable



### Reusable




### Extra filter




## PATIENT CASSETTE

### Reusable



### Patient cassette replacement parts



## Contents

- Starter kit consumables
- Heat and Moisture Exchange (HME)
- Patient connection
- Anesthesia mask
- Viral/bacterial filter
- Manual ventilation
- O<sub>2</sub> flowmeter / Suction equipment
- Gas analyze
- Gas evacuation
- CO<sub>2</sub> absorber
- Patient cassette

### For local contact:

Please visit our website  
[www.maquet.com](http://www.maquet.com)

**MAQUET**  
GETINGE GROUP

Maquet Critical Care AB  
SE-171 54 Solna  
SWEDEN  
Phone: +46 (0) 8 730 73 00  
[www.maquet.com](http://www.maquet.com)

Products may be pending regulatory approvals to be marketed in your country. Contact your MAQUET representative for more information.

GETINGE GROUP is a leading global provider of products and systems that contribute to quality enhancement and cost efficiency within healthcare and life sciences. We operate under the three brands of ArjoHuntleigh, GETINGE and MAQUET. ArjoHuntleigh focuses on patient mobility and wound management solutions. GETINGE provides solutions for infection control within healthcare and contamination prevention within life sciences. MAQUET specializes in solutions, therapies and products for surgical interventions, interventional cardiology and intensive care.