


120x60cm			
ITEM	QTY	TITLE	NUMBER
19	1	Bordplade	64112061
20	1	Overramme	20-60800-101
21	1	Frontprofil	20-60320-108
22	1	Aksel	20-60800-112
23	1	Bundlæg	20-60800-105

90x60cm			
ITEM	QTY	TITLE	NUMBER
19	1	Bordplade	64109061
20	1	Overramme	20-60800-117
21	1	Frontprofil	20-60320-077
22	1	Aksel	20-60800-111
23	1	Bundlæg	20-60800-095

20-60320-140	2	Vippebeslag/hængsel Vision			18
93000068	3	Motorkabel			17
93000179	1	Kontrollboks CBD6S			16
93000175	2	Ben med motor			15
20-60800-119	1	Bøsning			14
96100916	4	Stilleskrue M10x30			13
93000182	1	Betjeningspanel vinkelregulering			12
93000039	1	Betjeningspanelop/ned			11
93000006	2	Netkabel			10
98000-555	1	Stopring alu $\varnothing 15/\varnothing 8,1$ L=6, M4			9
20-60800-093	1	Beslag 1 til motor			8
93000180	1	Motor til vinkelregulering			7
20-60800-098	1	Beslag 2 til motor			6
20-60800-090	1	Kardanled med $\varnothing 3$ huller			5
20-60800-088	1	Kobling			4
20-60320-026	1	Skråstillsingsbeslag til Vision			3
96000656	1	Gearhus for skråstilling			2
20-60800-078	2	Fod til bord			1
PART NUMBER	QTY	TITLE	L	B	ITEM



High/low B 90/120cm Reservevedele		Tegn. Sj	15-09-2014
		DS/ISO 2768 m	
ROPOX Better ways to better days Ringstedgade 221, 4700 Næstved Tlf.: 55750500 Fax.: 55750550		Vægt:	N/A N/A
		Tegning nr.: 20-60801-2	