

BIOPSY FORCEPS

POLYPECTOMY SNARES



BIOPSY FORCEPS



Precise and safe tissue removal

- ✓ High precision cut
- ✓ Coated
- ✓ Incl. tissue sampler



Standard biopsy forceps

Precise and safe tissue removal



COATED – OVAL CUPS

Features

- ☐ Oval cups
- ☐ High precision cut avoids risk of tissue damage
- ☐ Position markings
- ☐ Integrated ergonomic handle
- ☐ Incl. tissue sampler
- ☐ Individually sterile packed



Outer diameter	Working length	Cup length	Cup depth	Cup volume	Pieces per box	Manufacturer Product Code	
						Without spike	With spike
PEDIATRIC ENDOSCOPES FROM 2 MM							
1.8 mm	120 cm	1.93 mm / 2.45 mm	0.68 mm	2.19 mm ³	10	BF11-ECO-18-120-SU	BF12-ECO-18-120-SU
1.8 mm	160 cm	1.93 mm / 2.45 mm	0.68 mm	2.19 mm ³	10	BF11-ECO-18-160-SU	BF12-ECO-18-160-SU
FOR BRONCHOSCOPES WITH WORKING CHANNELS FROM 2 MM							
1.8 mm	180 cm	1.93 mm / 2.45 mm	0.68 mm	2.19 mm ³	10	BF11-ECO-18-180-SU	BF12-ECO-18-180-SU
1.8 mm	230 cm	1.93 mm / 2.45 mm	0.68 mm	2.19 mm ³	10	BF11-ECO-18-230-SU	BF12-ECO-18-230-SU
FOR ENTEROSCOPES WITH WORKING CHANNELS FROM 2 MM							
1.8 mm	280 cm	1.93 mm / 2.45 mm	0.68 mm	2.19 mm ³	10	BF11-ECO-18-280-SU	BF12-ECO-18-280-SU
FOR GASTROSCOPES AND COLONOSCOPES WITH WORKING CHANNELS FROM 2.5 MM							
2.3 mm	160 cm	2.5 mm / 3 mm	0.85 mm	5.37 mm ³	10	BF11-ECO-23-160-SU	BF12-ECO-23-160-SU
2.3 mm	180 cm	2.5 mm / 3 mm	0.85 mm	5.37 mm ³	10	BF11-ECO-23-180-SU	BF12-ECO-23-180-SU
2.3 mm	230 cm	2.5 mm / 3 mm	0.85 mm	5.37 mm ³	10	BF11-ECO-23-230-SU	BF12-ECO-23-230-SU



Coated



ECO biopsy forceps

Precise and safe tissue removal



COATED – OVAL CUPS

Features

- ☐ Oval cups
- ☐ High precision cut avoids risk of tissue damage
- ☐ Position markings
- ☐ Integrated ergonomic handle
- ☐ Individually sterile packed



Outer diameter	Working length	Cup length	Cup depth	Cup volume	Pieces per box	Manufacturer Product Code	
						Without spike	With spike
FOR GASTROSCOPES AND COLONOSCOPES WITH WORKING CHANNELS FROM 2.5 MM							
2.3 mm	180 cm	2.5 mm / 3 mm	0.85 mm	5.37 mm ³	10	BF11-ECO-S-23-180-SU	BF12-ECO-S-23-180-SU
2.3 mm	230 cm	2.5 mm / 3 mm	0.85 mm	5.37 mm ³	10	BF11-ECO-S-23-230-SU	BF12-ECO-S-23-230-SU



Jumbo biopsy forceps

Precise and safe tissue removal



COATED – OVAL CUPS

Features

- ☐ Oval cups with double cup volume compared to standard forceps
- ☐ High precision cut avoids risk of tissue damage
- ☐ Position markings
- ☐ Integrated ergonomic handle
- ☐ Incl. tissue sampler
- ☐ Individually sterile packed



Outer diameter	Working length	Cup length	Cup depth	Cup volume	Pieces per box	Manufacturer Product Code	
						Without spike	With spike
FOR COLONOSCOPES WITH WORKING CHANNELS FROM 3.2 MM							
3 mm	230 cm	4.0 mm / 3.6 mm	1.1 mm	8.79 mm ³	10	BF11-ECO-30-230-SU	BF12-ECO-30-230-SU





Alligator biopsy forceps

Precise and safe tissue removal



COATED – ALLIGATOR CUPS

Features

- ☐ Alligator cups
- ☐ High precision cut avoids risk of tissue damage
- ☐ Position markings
- ☐ Integrated ergonomic handle
- ☐ Incl. tissue sampler
- ☐ Individually sterile packed



Outer diameter	Working length	Cup length	Cup depth	Cup volume	Pieces per box	Manufacturer Product Code	
						Without spike	With spike
FOR BRONCHOSCOPES WITH WORKING CHANNELS FROM 2 MM							
1.8 mm	180 cm	2.05 mm / 2.7 mm	0.83 mm	2.22 mm ³	10	BF11-ALI-18-180-SU	BF12-ALI-18-180-SU
1.8 mm	230 cm	2.05 mm / 2.7 mm	0.83 mm	2.22 mm ³	10	BF11-ALI-18-230-SU	BF12-ALI-18-230-SU
FOR ENTEROSCOPES WITH WORKING CHANNELS FROM 2 MM							
1.8 mm	280 cm	2.05 mm / 2.7 mm	0.83 mm	2.22 mm ³	10	BF11-ALI-18-280-SU	BF12-ALI-18-280-SU
FOR GASTROSCOPES AND COLONOSCOPES WITH WORKING CHANNELS FROM 2.5 MM							
2.3 mm	180 cm	2.88 mm / 3.08 mm	0.98 mm	4.85 mm ³	10	BF11-ALI-23-180-SU	BF12-ALI-23-180-SU
2.3 mm	230 cm	2.88 mm / 3.08 mm	0.98 mm	4.85 mm ³	10	BF11-ALI-23-230-SU	BF12-ALI-23-230-SU



POLYPECTOMY SNARES



Tissue-conserving removal of polyps

- ✓ Various shapes
- ✓ Multifilament
- ✓ HF ablation



Company GmbH



Multifilament

mdd

Polypectomy snares


Tissue-conserving removal of polyps



HF ABLATION

Features


- Easy capturing of polyps that are difficult to reach
- Controlled and continuously variable incision thanks to dimensionally stable and bending-resistant multifilament snare wire
- Integrated ergonomic handle for comfortable use
- Olympus HF connection
- Individually sterile packed



Catheter outer diameter	Working length	Open size	Cutting wire diameter	Pieces per box	Manufacturer Product Code		
					Sickle	Hexagonal	Oval
FOR GASTROSCOPES AND COLONOSCOPES WITH WORKING CHANNELS FROM 2.6 MM							
2.4 mm	160 cm	15 mm	0.4 mm	10	CRE-DS-23-15-160	HEX-DS-23-15-160	OVA-DS-23-15-160
2.4 mm	160 cm	20 mm	0.4 mm	10	CRE-DS-23-20-160	HEX-DS-23-20-160	OVA-DS-23-20-160
2.4 mm	160 cm	25 mm	0.4 mm	10	CRE-DS-23-25-160	HEX-DS-23-25-160	OVA-DS-23-25-160
2.4 mm	160 cm	30 mm	0.4 mm	10	CRE-DS-23-30-160	HEX-DS-23-30-160	OVA-DS-23-30-160
2.4 mm	230 cm	10 mm	0.4 mm	10	CRE-DS-23-10-230	HEX-DS-23-10-230	OVA-DS-23-10-230
2.4 mm	230 cm	15 mm	0.4 mm	10	CRE-DS-23-15-230	HEX-DS-23-15-230	OVA-DS-23-15-230
2.4 mm	230 cm	20 mm	0.4 mm	10	CRE-DS-23-20-230	HEX-DS-23-20-230	OVA-DS-23-20-230
2.4 mm	230 cm	25 mm	0.4 mm	10	CRE-DS-23-25-230	HEX-DS-23-25-230	OVA-DS-23-25-230
2.4 mm	230 cm	30 mm	0.4 mm	10	CRE-DS-23-30-230	HEX-DS-23-30-230	OVA-DS-23-30-230
2.4 mm	230 cm	35 mm	0.4 mm	10	CRE-DS-23-35-230	HEX-DS-23-35-230	OVA-DS-23-35-230
2.4 mm	230 cm	50 mm	0.4 mm	10	-	-	OVA-DS-23-50-230



ROTATABLE POLYPECTOMY SNARES



**Easier capturing of polyps thanks
to rotatable snare**

- ✓ 360° rotation
- ✓ Multifilament
- ✓ HF ablation



Company GmbH



Rotatable

**Rotatable polypectomy snares**

Easier capturing of polyps thanks to rotatable snare

**HF ABLATION****Features**

- ☑ 360° rotation makes it easier to capture polyps
- ☑ Gentle tissue removal
- ☑ Easy capture and removal even of polyps that are difficult to reach
- ☑ Controlled and continuously variable incision thanks to dimensionally stable and bending-resistant multifilament snare wire
- ☑ Integrated ergonomic handle for comfortable use
- ☑ Olympus HF connection
- ☑ Individually sterile packed



Catheter outer diameter	Working length	Open size	Cutting wire diameter	Pieces per box	Manufacturer Product Code
					Oval – rotatable
FOR GASTROSCOPES AND COLONOSCOPES WITH WORKING CHANNELS FROM 2.6 MM					
2.4 mm	160 cm	15 mm	0.4 mm.	10	OVA-DR-24-15-160
2.4 mm	160 cm	20 mm	0.4 mm.	10	OVA-DR-24-20-160
2.4 mm	160 cm	25 mm	0.4 mm.	10	OVA-DR-24-25-160
2.4 mm	160 cm	30 mm	0.4 mm.	10	OVA-DR-24-30-160
2.4 mm	230 cm	10 mm	0.4 mm.	10	OVA-DR-24-10-230
2.4 mm	230 cm	15 mm	0.4 mm.	10	OVA-DR-24-15-230
2.4 mm	230 cm	20 mm	0.4 mm.	10	OVA-DR-24-20-230
2.4 mm	230 cm	25 mm	0.4 mm.	10	OVA-DR-24-25-230
2.4 mm	230 cm	30 mm	0.4 mm.	10	OVA-DR-24-30-230
2.4 mm	230 cm	50 mm	0.4 mm.	10	OVA-DR-24-50-230



MULTI-STAGE POLYPECTOMY SNARES



Up to three different sizes in one snare

- ✓ Multi-stage and 360° rotation
- ✓ Multifilament
- ✓ Cold ablation



Company GmbH



Multi-stage polypectomy snares


Up to three different sizes in one snare



COLD ABLATION

Features

- ☐ The multi-stage design combines three sizes and shapes in one snare, which makes it unnecessary to replace the snare during a treatment
- ☐ The design with additional 360° makes it easier to capture polyps
- ☐ Gentle tissue removal
- ☐ Easy capture and removal even of polyps that are difficult to reach
- ☐ Controlled incision thanks to dimensionally stable and bending-resistant multifilament snare wire
- ☐ Integrated ergonomic handle for comfortable use
- ☐ Individually sterile packed



Outer catheter diameter	Working length	Open size	Wire thickness	Pieces per box	Manufacturer Product Code	
					Multi-stage	Multi-stage – rotatable
2.4 mm	230 cm	10 mm	0.25 mm	10	MST-CS-24-02-10-230	MST-CR-24-02-10-230
2.4 mm	230 cm	10 mm	0.30 mm	10	MST-CS-24-03-10-230	MST-CR-24-03-10-230



HYBRID POLYPECTOMY SNARES



Cold and HF ablation between 3 and 15 mm

- ✓ Multi-stage and 360° rotation
- ✓ Multifilament
- ✓ Cold and HF ablation




Multi-stage and rotatable

**Hybrid polypectomy snares**

Cold and HF ablation between 3 and 15 mm

**COLD AND HF ABLATION****Features**

- ☐ Various opening widths for ablating polyps between 3 and 15 mm with and without HF current
- ☐ 360° rotation makes it easier to capture polyps
- ☐ Gentle tissue removal
- ☐ Each positioning for flat lesions thanks to braided snare
- ☐ Controlled incision thanks to dimensionally stable and bending-resistant multifilament snare wire
- ☐ Integrated ergonomic handle for comfortable use
- ☐ Individually sterile packed



Outer catheter diameter	Working length	Open size	Wire thickness	Pieces per box
2.4 mm	230 cm	10 mm	0.4 mm	10
2.4 mm	230 cm	15 mm	0.4 mm	10

Manufacturer Product Code

Multi-stage – rotatable

OVA-HYB-24-10-230

OVA-HYB-24-15-230



Company GmbH



UpViser Aps
Adelgade 45, 1
8660 Skanderborg, Denmark
Tel.: +45 6171 7625
E-Mail: salg@upviser.com
www.upviser.dk



Exclusive Partner Europe
Body Products GmbH

mdd Company GmbH
Rengglochstr. 38
CH-6012 Obernau, Switzerland
Tel: +39 348 0026970
E-Mail: gb@mddcompany-gmbh.com
www.mddcompany-gmbh.com