



ENDOSCOPE TRANSPORT DRAWSTRING BAG

Easy transport of processed and used endoscopes

- ✓ Large volume
- ✓ Inspection window
- ✓ Moisture-proof fleece padding





Safe and protected transport



Endoscope transport drawstring bag


Easy transport of processed and used endoscopes



The endoscope drawstring bags were specially developed for transporting used endoscopes quickly and practically. The disposable drawstring bag meets the minimum requirements for transporting endoscopes safely from point A to point B in hospitals and doctor's offices.

Absorbent fleece padding

Any residual liquid from the used endoscope is safely collected.



Drawstring fastener

Two practical drawstrings for closing and carrying the bag easily



Inspection window

For a quick look at the contents in the drawstring bag



Saves space when folded up



Open it up and place the endoscope inside



To close, pull the drawstring

Features

- ☑ Large volume, making it also suitable for longer, flexible endoscopes
- ☑ Made out of polyethylene (PE), thickness 100 µm
- ☑ Two drawstrings for closing and carrying the bag easily
- ☑ Moisture-proof fleece padding
- ☑ Transparent blue inspection window
- ☑ Size when unfolded: 100 x 77 cm
- ☑ 100 % recyclable

Size L x W

Colour

Packaging unit

Manufacturer Product Code

100 x 77 cm

Blue

100

ETB-001-CL