

# Splentis

POP TISSUE ANCHORING SYSTEM

RELIABLE APICAL SUSPENSION



### Splentis POP Tissue Anchoring System

is a surgical kit consisting of polypropylene implants for SSL (Sacrospinous Ligament) fixation through vaginal approach. The kit is made up of TAS (Tissue Anchoring System) anchors and a Retractable Insertion Guide.

### Splentis POP Tissue Anchoring System

is indicated for Colpopexy with SSL fixation in minimal invasive vaginal surgery, through anterior or posterior approach.

### Simple and secure procedure



- Vaginal approach, anterior or posterior.
- Fixation points at SSL.
- Suspension points at pericervical ring, uterosacral ligament or vaginal vault.

## UNIQUE FIXATION SYSTEM: TAS ANCHORS

- High penetration force and accuracy.
- Reliable fixation.
- Safety in relation to vascular and neural structures.

### BENEFITS

- Certainty and simplicity for SSL fixation.
- Durable results.
- Low complication rate.



6 Fixation points at 360°

Safety stop

Reference:

1. Biosafety Assesment. Data on file.

## ACCURATE IMPLANTATION: RETRACTABLE INSERTION GUIDE - RIG

- Ergonomic design.
- Easy access to SSL through anterior or posterior approach.
- Retractable mechanism for connecting TAS anchors.

### BENEFITS

- Precision and safety in surgical maneuvers.
- Possibility to choose the right approach for each patient.
- Total control of retention and release of TAS anchors.




### Description

- 3 TAS (Tissue Anchoring System) anchors with sutures.
- 1 Retractable Insertion Guide.

The kit also includes:

- Tissue Reinforcement implant for suspension sites.
- 1 needle.



# Splentis

POP TISSUE ANCHORING SYSTEM

## ORDERING INFORMATION

SPLENTIS POP TISSUE ANCHORING SYSTEM

ORDER NUMBER: KIT –UT-01

- 3 TAS (Tissue Anchoring System) anchors with sutures: Ref: TAS
- 1 Retractable Insertion Guide: Ref: DPN-MNL
- 1 Reinforcing implant: Ref: MSAP
- 1 Eye suture needle: Ref.: ESN

B-90-22 (01) / 01-Mar-13



**Promedon**  
*People + Innovation*

[www.promedon.com](http://www.promedon.com)

CE 0197