

# DIGITAL URETEROSCOPE CONTROL UNIT



## Specifications and Technical Data

<b>Power Supply</b>	Voltage 100-240V AC Frequency 50/60 Hz
<b>CMOS Resolution</b>	400 (horizontal) x 400 (Vertical) pixels
<b>Outputs</b>	HDMI, DVI, USB3.0
<b>Examination Lamp</b>	2 LEDs
<b>Image Function</b>	white balance, brightness adjustment, zoom in & out, photo and video replay, foot panel for photo and freeze
<b>Image Storage</b>	32G SD card for photo and video

### Shanghai Anqing Medical Instrument Co.,Ltd.

Room 202, 2nd Floor, No.7 Building&Room 202, 2nd Floor, No.4 Building, 150Cailun Road, Zhangjiang High-Tech Park, 201210 Shanghai, China  
TEL: 8621 50191132 FAX: 8621 50191132-8004 ZIP: (021)201210  
Email:global@innovexmed.com.cn

# Single-Use Flexible Ureteroscope



A NEW SCOPE FOR EVERY PATIENT

EU-Scope™

MODEL:US31B-12

# No Contamination No Repair No Sterilization



## Product Features:

- 1 Field of view: 120°
- 2 Depth:3MM- 50MM
- 3 160K Resolution
- 4 Passive tip bending
- 5 Device locking function
- 6 Mutiple outputs: HDMI, DVI, USB 3.0



### Easy maneuverability

275° up & down deflection and precise tip motion makes navigation easy

### High quality image

2.5X High imaging pixel guarantees clear & sharp images during the lithotripsy

