


## SPECIAL PRODUCTS

### FLEXIBLE INJECTION NEEDLE

For injection of substances into the lower urinary tract

#### DESCRIPTION



#### Flexible Injection Needle (Injectra)

- ◆ Safe needle guidance and needlestick protection due to interlocking mechanism
- ◆ Ergonomic handle
- ◆ Medical stainless steel needle with 15° grinding insures optimal puncture characteristics


Size (FR/CH)	Needle Diameter (G)	Needle Length (mm)	Working length (cm)	Order Code
06	22	04	66	FIN-57-066

## SPECIAL PRODUCTS

### FLEXIBLE INJECTION NEEDLE

For injection of substances into the lower urinary tract

#### DESCRIPTION



#### Short flexible Needle

- ◆ Tuohy Borst Adapter for fixation on cystoscope working channel and for splash protection
- ◆ Easy and continuous orientation of needle position (grinding) via marking line

Size (FR/CH)	Needle Diameter (G)	Needle Length (mm)	Working length (cm)	Order Code
Flexible Needle (cylinder clearance 0,9 ml)				
05	23	08	38	FIN-05-380
06	21	08	38	FIN-06-380
Flexible Needle (cylinder clearance 0,9 ml) with Tuohy Borst Adapter				
05	23	08	36	FIN-TB-05-360
06	21	08	36	FIN-TB-06-360