

TABLE OF CONTENT



URETERAL STENTS
Chapter 1




NEPHROSTOMY
Chapter 2



SUPRAPUBIC BLADDER
DRAINAGE
Chapter 3



ACCESSORIES
Chapter 4



URETERAL CATHETERS
Chapter 5



FOLEY CATHETERS
Chapter 6



STONE BASKETS
Chapter 7



URODYNAMIC
Chapter 8




SPECIAL PRODUCTS
Chapter 9



PAEDIATRIC UROLOGY
Chapter 10


STONE BASKETS



STONE BASKETS

High elastic Nitinol material, with detachable handle acc. to Lahme

“TIPLESS”




- ♦ Delicate 3-, 4- or 6-wire-baskets without tip for the flexible uretero-rensoscopy
- ♦ For grasping of smallest stone concretments also on tissue edge
- ♦ With Tuohy Borst Adapter to protect the basket tip

Size (FR/CH)	Wires	Working Length (cm)	Basket Length (mm)	Basket Ø (mm)	Basket Shape	Order Code
1.2	3	120	12	8	straight	STN-T123-120
1.8	4	120	15	11	straight	STN-T184-120
1.8	4	120	15	11	helical	STN-HT184-120
1.8	6	120	13	9	straight	STN-T186-120
2.2	4	120	18	11	straight	STN-T224-120


Handle acc. to Lahme for all stone baskets

- ♦ Dis- and reconnectable handle: endoscope can be removed while basket is situated. Handle and basket can get attached afterwards
- ♦ Fixation of the stone in the basket due to spring mechanism
- ♦ Ergonomic handle design with non-slip surface
- ♦ Basket opens and closes within image area: no adjustment of endoscope necessary



STONE BASKETS

“HELICAL”




Size (FR/CH)	Wires	Working Length (cm)	Basket Length (mm)	Basket Ø (mm)	Tip length (mm)	Order Code
2.5	4	90	30	12	3	STN-L254-090
2.5	4	120	30	12	3	STN-L254-120
3.0	4	90	31	12	3	STN-L304-090
3.0	4	120	31	12	3	STN-L304-120
4.0	4	90	37	15	3	STN-L404-090

“FILIFORM”



Size (FR/CH)	Wires	Working Length (cm)	Basket Length (mm)	Basket Ø (mm)	Tip length (mm)	Order Code
3.0	4	90	30	12	40	STN-F304-090

“STRAIGHT”



Size (FR/CH)	Wires	Working Length (cm)	Basket Length (mm)	Basket Ø (mm)	Tip length (mm)	Order Code
2.5	4	120	30	12	3	STN-G254-120
3.0	4	120	31	12	3	STN-G304-120